

# Macawi Respiratory Module pneumatic evaluation kit.

RS232 communication cable





Pneumatic safety valve



Power cable (24V)



Macawi Respiratory Module (100)



Patient flow sensor  
(adult, ped/neo on request)



Corner piece  
connector



Pneumatic expiration valve



Expiration side  
plug



Patient side  
plug



Tube connector



T-piece connector



USB drive  
(IFU / documentation)

## EVALUATION KIT OF THE MACAWI RESPIRATORY MODULE

To give you a headstart in your ventilator development, we have created an evaluation kit of the Macawi Respiratory Module (MRM). With this plug & play evaluation kit you can easily test the performance of the respiratory module and adapt relevant ventilation parameters during the indesign engineering process.

With the evaluation kit, the following will be available for you:

- Access to the product support page on our website
- A technical kick-off session with our customer support engineer
- Continued online technical support during your development process regarding the MRM

To start testing with the evaluation kit, you need a computer with Windows OS, hose set, test lung and a 24V DC power supply. These items are not included in the evaluation kit.